

# RITMSCENAR Professional Device Specifications

RITMSCENAR® Super Pro

RITMSCENAR® Life Master

RITMSCENAR® Pro Plus

RITMSCENAR® Pro

FUNCTION / DEVICE



Display	Color LCD display	B & W LCD Display	B & W LCD Display	B & W LCD Display
Treatment Timer	Off, 5m, 1 0m, 15m, 20m, 30m	N/A	N/A	N/A
Damping	Off,1,2,3,4,VAR	Sc 1, 2, 3, 4, VAR	Sc 1, 2, 3, 4, VAR	Sc 1, 2, 3, 4, VAR
Amplitude Modulation	"1 :1", "2:1 ", "3:1 ", "4:1 ", "5:1", BEE	"1 :1", "2:1 ", "3:1 ", "4:1 ", "5:1", BEE	"1 : 1", "2 : 1", "3 : 1", "4 : 1", "5 : 1"	"1 : 1", "2 : 1", "3 : 1", "4 : 1", "5 : 1"
Frequency Modulation	From 30 to 120 Hz, with (7±2) sec cycle	From 30Hz to 120Hz, with (7±2) sec cycle	From 30Hz to 120Hz, with (7±2) sec cycle	From 30Hz to 120Hz, with (7±2) sec cycle
Swing Mode	Standard SW1, SW2, SW3, SW4	Standard SW1, SW2, SW3, SW4	Standard SW1, SW2, SW3, SW4	
Bio Gap	Yes	N/A	N/A	N/A
Bio Frequency	Yes	N/A	N/A	N/A
Influencing Frequencies	From 0,674Hz to 500Hz	From 0,641Hz to 350Hz	From 15Hz to 350Hz	From 15Hz to 350Hz
Intensity	From 1 to 8	From 1 to 8	From 1 to 8	From 1 to 8
Audio Schemes	1	1	1	1
Dr. <b>Voll's</b> frequencies	<b>128 Voll's frequencies (from 0.641 up to 14.9 Hz)</b> were added to usual-al frequencies.	<b>128 Voll's frequencies (from 0.641 up to 14.9 Hz)</b> were added to usual-al frequencies.	N/A	N/A
Dosing regimens:	"Diag1" integrally dose and zero. "Diag 2" is a differential dose	"Diag1" integrally dose and zero. "Diag 2" is a differential dose	"Diag1" integrally dose and zero. "Diag 2" is a differential dose	"Diag1" integrally dose and zero. "Diag 2" is a differential dose
Screening Mode: D-5	3 modes: "Diag3", "Diag4", "Diag5" Dose 5 is used to carry out electro puncture Diagnostics in 24 points (6 points on each extremity) with a local electrode.	3 modes: "Diag3", "Diag4", "Diag5" Dose 5 is used to carry out electro puncture Diagnostics in 24 points (6 points on each extremity) with a local electrode.	N/A	N/A
Additional modes:	6 additional ranges	N/A	N/A	N/A
Screen Backlight Timer	Off, 1 sec, 2 sec, 3 sec, 5 sec, 10 sec, 30 sec	Off, 1 sec, 2 sec, 3 sec, 5 sec, 10 sec, 30 sec	Off, 1 sec, 2 sec, 3 sec, 5 sec, 10 sec, 30 sec	Off, 1 sec, 2 sec, 3 sec, 5 sec, 10 sec, 30 sec
Standby Mode	Off, 5min, 10min, 15min, 20min, 30min	On, 5min, 10min, 20min, 5min, 10min, 20min	On, 5min, 10min, 20min, 5min	On, 5min, 10min, 20min, 5min
Volume Settings	Mute 1, 2, 3	Mute 1, 2, 3	Mute 1, 2, 3	Mute 1, 2, 3
Warranty	2 years	2 years	2 years	2 years